

Specification	Imperial	Metric
Safe Working Load	341 lbs	155 kgs
Maximum Overall Length	51.0 inches	1300 mm
Minimum Overall Length	49.0 inches	1250 mm
Maximum Overall Height	73.0 inches	1860 mm
Minimum Overall Height	53.5 inches	1360 mm
Height (Folded)	17.1 inches	450 mm
Depth (Folded)	46.4 inches	1180 mm
Width	21.6 inches	550 mm
Spreader Bar - Maximum Height	66.5 inches	1690 mm
Spreader Bar - Minimum Height (usable)	15.3 inches	390 mm
Height at Maximum Reach	46.0 inches	1170 mm
Reach at Maximum Height	25.6 inches	650 mm
Reach at Minimum Height	13.7 inches	350 mm
Maximum Reach (centre of spreader bar to front of the mast)	32.8 inches	835 mm
Turning Radius	55.9 inches	1420 mm
Legs Open - External Width	42.5 inches	1080 mm
Legs Open - Internal Width	39.3 inches	1000 mm
Legs Closed - External Width	26.3 inches	670 mm
Legs Closed - Internal Width	22.4 inches	570 mm
Overall Height of Legs	4.5 inches	115 mm
Ground Clearance	1.1 inches	30 mm
Front Twin Castors	3.0 inches	75 mm
Rear Braked Castors	3.9 inches	100 mm

Weights	Imperial	Metric
Mast, Base & Boom Assembly	63.1 lbs	28.9 kgs
Power Pack	6.6 lbs	3.0 kgs
Total	69.7 lbs	31.7 kgs
Base Assembly	32.8 lbs	14.9 kgs
Mast & Boom Assembly (not including battery)	30.9 lbs	14.0 kgs

Specifications may vary due to manufacturing tolerances

#### Electrical Specifications

Battery 1 X 12 volt rechargeable sealed lead acid type  
 Battery capacity 3.2 Ampere hours  
 Charger rated input 230Vac 50/60Hz  
 Charger rated output 27.4/29.0 VDC@0.8A

#### Warranty

The Advance comes with a 2 year warranty as standard.  
 Please refer to your user manual for further information.

#### Electrical Shock Protection

Charger - class II  
 Lift - internal power source

Degree of shock protection

Charger - Type B

Lift - Type B

Intended operating environment : >+5° <+40°

Outside this environment functionality and safety may be compromised



# Oxford®



Joerns Healthcare Limited  
 Drakes Broughton Business Park, Worcester Road, Drakes Broughton  
 Pershore, Worcestershire, WR10 2AG United Kingdom  
 (T) +44 (0) 844 811 1156 (F) +44 (0) 844 811 1157  
 www.joerns.co.uk • email: info@joerns.co.uk