

Specification	Imperial	Metric
Safe Working Load	440 lbs	200 kgs
Maximum Overall Length	43.3 inches	1100 mm
Minimum Overall Length	41.0 inches	1040 mm
Maximum Overall Height	66.9 inches	1700 mm
Minimum Overall Height	48.4 inches	1230 mm
Maximum Height to Attachment Point	64.6 inches	1640 mm
Minimum Height to Attachment Point	37.2 inches	945 mm
Turning Radius	49.2 inches	1250 mm
Legs Open - External Width	39.1 inches	995 mm
Legs Open - Internal Width	34.3 inches	870 mm
Legs Closed - External Width	24.8 inches	630 mm
Legs Closed - Internal Width	20.0 inches	510 mm
Widest Point (base width)	25.2 inches	640 mm
Overall Height of Legs	4.7 inches	120 mm
Ground Clearance	1.4 inches	35 mm
Front Twin Castors	4.0 inches	100 mm
Rear Braked Castors	4.0 inches	100 mm

Weights	Imperial	Metric
Mast, Base & Boom Assembly	130.1 lbs	59.0 kgs
Power Pack	6.6 lbs	3.0 kgs
Total	136.7 lbs	62.0 kgs
Base Assembly (not including battery)	40.8 lbs	18.5 kgs
Mast & Boom Assembly (not including battery)	66.1 lbs	30.0 kgs
Foot Tray	23.1 lbs	10.5 kgs

Specifications may vary due to manufacturing tolerances

#### Electrical Specifications

Battery 24 volt rechargeable sealed lead acid type  
 Battery capacity 3.2 Ampere hours  
 Charger rated input 100 - 240V AC 50/60Hz  
 Charger rated output 29.5V DC, Max 19W

#### Warranty

The Elevate comes with a 2 year warranty as standard.  
 Please refer to your user manual for further information.

#### Electrical Shock Protection

Charger - Class II  
 Lift - Internal power source  
  
 Degree of shock protection  
 Charger - Type B  
 Lift - Type B  
 Intended operating environment: >+5° <+40°  
 Outside this environment functionality and safety may be compromised





Joerns Healthcare Limited  
 Drakes Broughton Business Park, Worcester Road, Drakes Broughton  
 Pershore, Worcestershire, WR10 2AG United Kingdom  
 (T) +44 (0) 844 811 1156 (F) +44 (0) 844 811 1157  
 www.joerns.co.uk • email: info@joerns.co.uk