

Specification	Imperial	Metric
Safe Working Load	341 lbs	155 kgs
Maximum Overall Length	39.4 inches	1000 mm
Minimum Overall Length	36.8 inches	935 mm
Maximum Overall Height	63.0 inches	1600 mm
Minimum Overall Height	44.5 inches	1130 mm
Height (Folded)	20.5 inches	525 mm
Depth (Folded)	41.0 inches	1040 mm
Widest Point (Between support handles)	26.8 inches	680 mm
Maximum Height to Attachment Point	60.2 inches	1530 mm
Minimum Height to Attachment Point	28.4 inches	720 mm
Turning Radius	45.1 inches	1145 mm
Legs Open - External Width	38.2 inches	970 mm
Legs Open - Internal Width	34.5 inches	875 mm
Legs Closed - External Width	27.0 inches	685 mm
Legs Closed - Internal Width	22.2 inches	590 mm
Overall Height of Legs	4.3 inches	110 mm
Ground Clearance	1.2 inches	30 mm
Front Twin Castors	3.0 inches	75 mm
Rear Braked Castors	4.0 inches	100 mm

Weights	Imperial	Metric
Mast, Base & Boom Assembly	82.4 lbs	37.4 kgs
Power Pack	6.2 lbs	2.8 kgs
Total	88.6 lbs	40.2 kgs
Base Assembly	30.2 lbs	18.2 kgs
Mast & Boom Assembly (not including battery)	42.3 lbs	19.2 kgs

Specifications may vary due to manufacturing tolerances.
Patent Pending.

Electrical Specifications

Battery 1 X 24 volt rechargeable sealed lead acid type
Battery capacity 3.2 Ampere hours
Charger rated input 100 - 240 V AC/ 50/60Hz
Charger rated output 29.5 V DC Max 19 W

Warranty

The Journey comes with a 2 year warranty as standard.
Please refer to your user manual for further information.

Electrical Shock Protection

Charger - class II
Lift - internal power source

Degree of shock protection

Charger - Type B

Lift - Type B

Intended operating environment : >+5° <+40°

Outside this environment functionality and safety may be compromised



Oxford®



Joerns Healthcare Limited
Drakes Broughton Business Park, Worcester Road, Drakes Broughton
Persnore, Worcestershire, WR10 2AG United Kingdom
(T) +44 (0) 844 811 1156 (F) +44 (0) 844 811 1157
www.joerns.co.uk • email: info@joerns.co.uk