

Specification	Imperial	Metric
Safe Working Load	500 lbs	227 kgs
Maximum Overall Length	54.7 inches	1390 mm
Minimum Overall Length	50.8 inches	1290 mm
Maximum Overall Height	83.8 inches	2130 mm
Minimum Overall Height	76.6 inches	1945 mm
Maximum Height to Attachment Point - 4 Point Cradle	61.0 inches	1550 mm
Minimum Height to Attachment Point - 4 Point Cradle	15.4 inches	390 mm
Maximum Height to Attachment Point - 6 Point Spreader Bar	77.7 inches	1975 mm
Minimum Height to Attachment Point - 6 Point Spreader Bar	26.6 inches	675 mm
Turning Radius	61.0 inches	1550 mm
Legs Open - External Width	53.5 inches	1360 mm
Legs Open - Internal Width	48.4 inches	1230 mm
Legs Closed - External Width	29.9 inches	760 mm
Legs Closed - Internal Width	25.1 inches	640 mm
Overall Height of Legs	4.3 inches	110 mm
Ground Clearance	1.2 inches	30 mm
Front Twin Castors	4.0 inches	100 mm
Rear Braked Castors	4.0 inches	100 mm

Weights	Imperial	Metric
Mast, Base & Boom Assembly	142 lbs	64.5 kgs
4 Point Positioning Cradle (manual operation)	13.2 lbs	6.0 kgs
Power Pack	6.6 lbs	3.0 kgs
Total	162 lbs	73.5 kgs
4 Point Positioning Cradle (powered operation)	19.8 lbs	9.0 kgs

Specifications may vary due to manufacturing tolerances

Electrical Specifications (Lift)

Battery 1 X 12 volt rechargeable sealed lead acid type
 Battery capacity 3.2 Ampere hours
 Charger rated input 230Vac 50/60Hz
 Charger rated output 27.4/29.0 VDC@0.8A

Electrical Specifications (Powered cradle only)

Battery 24V DC rechargeable NiMH sealed
 Battery capacity 1400 mAh
 Charger rated input 100 - 240V AC / 24 V DC 50/60Hz Max. 75mA
 Charger rated output Max 28.4V DC

Warranty

The Stature comes with a 2 year warranty as standard.
 Please refer to your user manual for further information.


Electrical Shock Protection

Charger - class II
 Lift - internal power source

Degree of shock protection
 Charger - Type B
 Lift - Type B
 Intended operating environment : >+5° <+40°
 Outside this environment functionality and safety may be compromised



Oxford®



Joerns Healthcare Limited
 High Street, Wollaston, Stourbridge, West Midlands DY8 4PS England
 (T) +44 (0)1384 446622 (F) +44 (0)1384 446601
 www.joerns.co.uk • email: info@joerns.co.uk